

Instrumentation Kits

BM Kit 1



BM Kit 2



Custom Kit



NIBP Service Kits



Fundamental Field Service Kit

This **Field Service Kit** was designed to adequately equip your field service representative with the fundamental testers, analyzers and simulators he will most likely need when visiting the average clinic. This **specially priced Kit** contains a Defibrillator Tester, Pressure Meter, ECG simulator and Medical Equipment Safety Analyzer. [More..](#)

Advanced Field Service Kit

This **Field Service Kit** was designed to adequately equip your field service representative with more advanced testers, analyzers and simulators he will most likely encounter when visiting Ambulatory Surgery Centers and other more complex off site areas. This **specially priced Kit** contains a Defibrillator Tester, Pressure Meter, two ECG simulators, external pacemaker tester and Medical Equipment Safety Analyzer. [More..](#)

Design the Kit you need

Create the **Field Service Kit** which is best suited to your own application. Your selected items will come with a durable carrying case and a **special price**. [More..](#)

SC-4 + OX Kit / SC-5 + OX Kit

SC-4 Kit - buy this kit and receive the complete SC-4 plus the battery boost option which eliminates the need for an AC adapter plus other extras. [More..](#)

SC-5 Kit - buy this kit and receive the complete SC-5 plus the battery boost option which eliminates the need for an AC adapter plus other extras. [More..](#)

FIELD SERVICE KIT

KIT-1 FUNDAMENTAL



Fundamental Field Service Kit

Basic needs for a Field Service Representative; includes defibrillator analyzer, safety analyzer, ECG simulator and pressure meter.

BM KIT-1

This **Field Service Kit** was designed to adequately equip your field service representative with the fundamental testers, analyzers and simulators he will most likely need when visiting the average clinic. This **specially priced Kit** contains a Defibrillator Tester, Pressure Meter, ECG simulator and Medical Equipment Safety Analyzer.

This Kit contains:

- ◆ **Delta 3000A** Tests all functions associated with defibrillation including AED, Bi-Phasic and Transthoracic Pacing
- ◆ **LKG 610** This safety analyzer conforms to both AAMI/ESI 1993, performs isolation testing and lead wire leakage, current, voltage, polarity, ground resistance and chassis leakage.
- ◆ **MiniSim 0.1** A basic 10 lead ECG simulator with normal sinus rhythm and arrhythmia waveforms.
- ◆ **DigiMano 1000** The 2000PSI is a pressure meter that can measure to 1999 mmHG or 19.99 PSI. The input uses the luer lok for fast connection.
- ◆ **Aluminum Travel Case with diced padding**



110 Toledo Street, Farmingdale, NY 11735, Phone: 631-531-0100, Fax: 631-531-0101

www.Netech.org **1-800-547-6557 (US)**

FIELD SERVICE KIT

KIT-2 ADVANCED



Advanced Field Service Kit

Created to meet the needs of a full functional senior Field Service Representative; includes defibrillator analyzer, safety analyzer, two ECG simulators and pressure meter.

BM KIT-2

This **Field Service Kit** was designed to adequately equip your field service representative with more advanced testers, analyzers and simulators he will most likely encounter when visiting Ambulatory Surgery Centers and other more complex off site areas. This **specialty priced Kit** contains a Defibrillator Tester, Pressure Meter, two ECG simulators, external pacemaker tester and Medical Equipment Safety Analyzer.

This Kit contains:

- ◆ **Delta 3000** Tests all functions associated with defibrillation including AED, Bi-Phasic and Transthoracic Pacing, as well as external pacemakers including AV sequential.
- ◆ **MultiPro 2000** This safety analyzer conforms to both AAMI/ESI 1993, performs isolation testing and lead wire leakage, current, voltage, polarity, ground resistance and chassis leakage. Also contains a full functional ECG Simulator.
- ◆ **MiniSim 1.3** Advanced 10 lead ECG simulator with normal sinus rhythm, arrhythmia waveforms, IBP waveforms, Temperature and respiration.
- ◆ **UniMano 800** The model 5000 is a pressure/vacuum meter that can measure from an absolute vacuum to 100 PSI to meet any challenge in the field. The input uses the luer lok for fast connection.
- ◆ Aluminum Travel Case with diced padding



110 Toledo Street, Farmingdale, NY 11735, Phone: 631-531-0100, Fax: 631-531-0101

www.Netech.org **1-800-547-6557 (US)**

FIELD SERVICE KIT

KIT-Custom



Custom Field Service Kit

Create your own test case from our wide selection of medical test equipment. Why leave such an important decision to someone else?

BM KIT-Custom

This **specially priced Kit** exists only in your mind. You decide what your needs are and we will make sure they work to the excellence you expect.

Take your pick from our wide range of test devices. Here is a sample of some of the equipment that can be selected:

Defibrillator Analyzers

Delta 1500 AED Tester
Delta 3000A-Defib Tester
Delta 3000 Defib/Pacer Tester

Safety Analyzers

LKG 601
LKG 610
MultiPro 2000

Simulators

MircoSim - Arrh
MircoSim - Pressure
MiniSim Basic 1000
MiniSim 1000
Holter Simulator
EEG Simulator
Cardiac Output

Pressure Meters

DigiMano 1000
DigiMano 2000
UniMano 800
UniMano 880

Flow Meter

Pacemaker Analyzers

PMT 100
PMT 2000

Call Your Sales Representative [1-631-531-0100](tel:1-631-531-0100) or contact Sales@Netech.ORG for assistance.



110 Toledo Street, Farmingdale, NY 11735, Phone: 631-531-0100, Fax: 631-531-0101

www.Netech.org **1-800-547-6557 (US)**